



ANALYTICAL BALANCE BFU1CD8

ANALYTICAL BALANCE BFU1CD8

TOUCH SCREEN BALANCE



1. Large touch screen
2. Powerful functions like weighing, parts counting, percent weighing, check weighing, animal weighing, density weighing, etc.
3. Standard anti-static module
4. 0.01mg maximum weighing range 100g
5. Electronic level adjustment indicator module
6. GLP/GMP

SPECIFICATIONS

Model	BFU1CD8
Weighing Range	0-50g/0-200g
Calibration	External
Readability	0.01/0.1mg
Stabilization Time	0.01mg≤13S, 0.1mg≤4S
Pan Size	Φ80mm
Dimension/ N.W.	360*220*350mm/8KG
Shipping Size/ G.W.	520*380*500mm/10kg
Alt Name	Touch Screen Balance



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com