





## ANALYTICAL BALANCE BFU1CG5

# ANALYTICAL BALANCE BFU1CG5

BALANCE



## SPECIFICATIONS

|                      |                                                                            |
|----------------------|----------------------------------------------------------------------------|
| Model                | BFU1CG5                                                                    |
| Weighing Range       | 0-20g/0-100g                                                               |
| Calibration          | External                                                                   |
| Readability          | 0.01/0.1mg                                                                 |
| Stabilization Time   | 0.01mg≤13S, 0.1mg≤4S                                                       |
| Pan Size             | Φ85mm                                                                      |
| Dimension/ N.W.      | 350*215*340mm/8KG                                                          |
| Shipping Size/ G.W.  | 520*380*500mm/10kg                                                         |
| Standard Accessories | Calibration Weight (for external calibration), Adapter, Manual, Dust cover |
| Alt Name             | 60 0.01mg Balance                                                          |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)