



## ANALYTICAL BALANCE BFU1CH20

# ANALYTICAL BALANCE BFU1CH20

TOUCH SCREEN BALANCE



1. 4.3-inch color touch screen
2. Powerful functions like weighing, parts counting, percent weighing, check weighing, animal weighing, density weighing, etc.
3. Electromagnetic balance sensor.
4. External/ Internal calibration.
5. Supports PC connection. Wireless bluetooth printing is optional.

## SPECIFICATIONS

|                    |                      |
|--------------------|----------------------|
| Model              | BFU1CH20             |
| Maximum Capacity   | 500g                 |
| Calibration        | Internal             |
| Readability        | 1mg                  |
| Stabilization Time | ≤4S                  |
| Pan Size           | Φ85mm                |
| Dimension/N.W.     | 360*220*350mm / 8kg  |
| Shipping Size/G.W. | 520*380*500mm / 10kg |
| Alt Name           | Touch Screen Balance |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)