



ANALYTICAL BALANCE BFU1CH3

ANALYTICAL BALANCE BFU1CH3

TOUCH SCREEN BALANCE



1. 4.3-inch color touch screen
2. Powerful functions like weighing, parts counting, percent weighing, check weighing, animal weighing, density weighing, etc.
3. Electromagnetic balance sensor.
4. External/ Internal calibration.
5. Supports PC connection. Wireless bluetooth printing is optional.

SPECIFICATIONS

Model	BFU1CH3
Maximum Capacity	200g
Calibration	External
Readability	0.1mg
Stabilization Time	≤4S
Pan Size	Φ85mm
Dimension/N.W.	360*220*350mm / 8kg
Shipping Size/G.W.	520*380*500mm / 10kg
Alt Name	Touch Screen Balance



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com