



COD ANALYZER BFU1T1

COD ANALYZER BFU1T1

COD ANALYZER



1. 7 inch high-resolution touch screen with friendly GUI.
2. The unique high-precision filter system makes the instrument more stable and accurate.
3. No blank calibration, automatic measurement, direct reading of concentration.
4. With built-in standard curves. Users can set up curves according to different measurement requirements.
5. Support zero point and up to 7-point calibration. Support calibration restoring function to restore factory calibration data.
6. With a large-capacity storage function, can store 5000 sets of test data.
7. Built-in thermal printers can print current data and stored data.

SPECIFICATIONS

| | |
|----------------------|---|
| Model | BFU1T1 |
| Measurement Method | Fast digestion spectrophotometry |
| Measuring Range | 10-10000mg/L |
| Accuracy | ≤5% |
| Resolution | 0.01 |
| Light Source | Imported Silicon Photodiode |
| Display | 3-inch color screen |
| Output | Micro USB |
| Power | Power adapter/Rechargeable lithium battery |
| Waterproof Level | IP54 |
| Ambient Temperature | (5-40) °C |
| Ambient Humidity | Relative humidity ≤ 85% (No condensation) |
| Dimensions/N.W. | 100x210x70 mm/1kg |
| Shipping size/G.W. | 355x265x150mm/2kg |
| Standard Accessories | Sample bottle x 2 pcs, Power adapter(Charger), Manual |
| Alt Name | COD Analyzer |



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com