



COLORIMETER BFU1V1

COLORIMETER BFU1V1

COLORIMETER



1. 7 inch high-resolution touch screen with friendly GUI.
2. The unique high-precision filter system makes the instrument more stable and accurate.
3. Direct reading of concentration.
4. With built-in standard curves. Users can set up curves according to different measurement requirements.
5. Support up to 7-point calibration. Support calibration restoring function to restore factory calibration data.
7. With a large-capacity storage function, can store 5000 sets of test data.
8. Built-in thermal printers can print current data and stored data.

SPECIFICATIONS

Model	BFU1V1
Measurement Method	Platinum-Cobalt Colorimetry
Measuring Range	0-500PCU
Abs Range	-0.3~3A
Colorimetric Method	Tube Colorimetric
Accuracy	≤5%
Resolution	0.01 PCU
Light Source	Imported Silicon Photodiode
Display	7-inch high-resolution 1024 * 600 touch screen
Output	USB (for data transfer)
Power	AC 220V/50Hz or AC 110V/60Hz
Ambient Temperature	(5-40) °C
Ambient Humidity	Relative humidity ≤ 85% (No condensation)
Dimensions/N.W.	460x380x200 mm/8kg
Shipping size/G.W.	520x450x300mm/10kg
Standard Accessories	10mm glass cuvette x 4 pcs, 30mm glass cuvette x 2 pcs, Power cord, Manual
Alt Name	Colorimeter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com