



CONDUCTIVITY METER BFU1Z2

CONDUCTIVITY METER BFU1Z2

PEN TYPE PH/CONDUCTIVITY METER



1. Mini design with size 180x50x20mm.
2. Auto calibration, auto power-off, low voltage alarm, electrode failure warning.
3. IP57 waterproof, replaceable electrode.

SPECIFICATIONS

Model	BFU1Z2
Description	pH Meter
Measuring Range	-1.00~15.00pH
Resolution	0.01pH
Accuracy	±0.05pH
Temp. Compensation	0 ~ 60°C
Calibration	1-3 Point
Power	Lithium batteries (CR2032) x 2
Humidity	0.85
Standard Accessories	Backup batteries
Alt Name	Pen Type pH/Conductivity Meter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com