

## MEDICAL AUTOCLAVE BGA1H1

# MEDICAL AUTOCLAVE BGA1H1

## CLASS N STEAM STERILIZER

Designed for sterilization of solid instruments, medical class N autoclave is a non vacuum sterilizer which removes air from the chamber by thermodynamic air displacement. Safety features with a door locking system, these autoclaves are ideal for small spaces.

## SPECIFICATIONS

Model	BGA1H1
Capacity	18 L
Product Size	450X550X350 mm
Chamber Dimension	φ249x350 mm
Package Size	670x550x500 mm
Gross Weight	30 Kg
Net Weight	23 Kg
Power	800 W
Voltage	AC220/110V+10%, 50/60 Hz+ 1Hz
Alt Name	Class N Steam Sterilizer

## APPLICATIONS

Medical, Dental, Laboratory, Pharmaceutical, Research.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)