



DRY BATH INCUBATOR BGH1B1 (BDIB-117)

DRY BATH INCUBATOR BGH1B1

DRY BATH INCUBATOR

Round metal blocks and more compact structure, suitable for most centrifuge tubes. Block replacement is convenient, easy to clean and disinfect.



SPECIFICATIONS

Model	BGH1B1
Old Model	BDIB-117
Temperature control range	RT. +5 °C-100 °C
Timer range	1 min-99 h 59 min
Temp. control accuracy	±0.5 °C
Heating time (25 °C to 100 °C)	≤15 min (@20 °C~30 °C)
Certificate	CE
Dimension (LxWxH) mm	225x160x197 mm
Weight	3 kg
Alt Name	Dry Bath Incubator

ACCESSORIES FOR PURCHASE

No	Name	Type	Description
1	Block	A-30	96x0.2 mL tube
2	Block	B-30	29x0.5 mL+ 32x1.5 mL tube
3	Block	C-30	61x0.5 mL tube
4	Block	D-30	41x1.5 mL tube
5	Block	E-30	41x2.0 mL tube
6	Block	F-30	61xΦ6 mm tube
7	Block	H-30	28xΦ12 mm tube
8	Block	I-30	28xΦ13 mm tube
9	Block	J-30	21xΦ15 mm tube
10	Block	K-30	21xΦ16 mm tube
11	Block	L-30	17xΦ18 mm tube
12	Block	M-30	14xΦ20 mm tube
13	Block	N-30	9xΦ25 mm tube
14	Block	O-30	4xΦ32 mm tube



1



2



3



4



5



6



7



8



9



10



11



12



13



14



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com