



MINI CENTRIFUGE BGH1O1 (BCMI-510)

MINI CENTRIFUGE BGH101

MINI CENTRIFUGE



Composite shell, shatter-resistant.
The electronic braking system for rapid sample processing.
Quick release rotor tub, allows for easy rotor exchange.
Two speed options: 4,000 & 6,000 rpm.
It includes 2 quick release interchangeable rotors.
--8x1.5 / 2.0 mL rotor --4x8x0.2 mL PCR strip rotor
Two speed options
Quick Release Rotor

SPECIFICATIONS

Model	BGH101
Old Model	BCMI-510
Rotor	8x1.5 / 2.0 mL rotor 4x8x0.2 mL PCR strip rotor 8 or 32x0.2 mL rotor 8x0.5 mL rotor
Speed	4000 rpm & 6000 rpm
Max. RCF	1,200xg (4000 rpm) 2,000xg (6000 rpm)
Power supply	AC 100-240 V, 50/60 Hz
Alt Name	Mini Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Description
1	Rotor	8x1.5 mL / 2.0 mL rotor
2	Strip rotor	4x8x0.2 mL PCR strip rotor
3	tube	Adapter for 0.5 mL tube
4	tube	Adapter for 0.2 mL tube

FEATURES

Composite shell, shatter-resistant.
The electronic braking system for rapid sample processing.
Quick release rotor tub, allows for easy rotor exchange.
Two speed options: 4,000 & 6,000 rpm.
It includes 2 quick release interchangeable rotors.
--8x1.5 / 2.0 mL rotor --4x8x0.2 mL PCR strip rotor



Two speed options



Quick Release Rotor



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com