



BIOLOGICAL INDICATOR INCUBATOR BGH1P2 (BIMA-112)

BIOLOGICAL INDICATOR INCUBATOR BGH1P2

BIOLOGICAL INDICATOR INCUBATOR

Biological indicator incubators are designed for temperature incubation of the EO biological indicators. Both of them are compact, tabletop incubators.



LCD display.

Users can set the incubating temperature and time.

Exact correspondence of the set and actual temperature.

It is easy to set up and use.

The Incubator has many wells for placement of the biological indicator vials and the wells can be custom-made.

SPECIFICATIONS

Model	BGH1P2
Old Model	BIMA-112
Temp. control range	RT. +5 °C~100 °C
Timer range	1 min~99 h 59 min
Temp. control accuracy	±0.5 °C
Display accuracy	±0.1 °C
Heating time (25 to 57 °C)	≤6 min
Dimension (DxWxH) mm	120x152x112 mm
Weight	0.85 kg
Temp. stability @57 °C	±0.35 °C
Alt Name	Biological Indicator Incubator

ACCESSORIES FOR PURCHASE

No	Name	Type	Description
1	Block	A	15xΦ10 mm
2	Block	B	15xΦ12 mm
3	Block	C	6xΦ15 mm
4	Block	D	6xΦ16 mm
5	Block	E	24xΦ8.8 mm
6	Block	F	18xΦ9 mm



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com