



FUNGAL GROWTH INCUBATOR BGH1Q1 (BIFG-501)

FUNGAL GROWTH INCUBATOR BGH1Q1

GEL CARD INCUBATOR

A gel card incubator can be suitable for many types of gel cards from different manufacturers.



2 independent incubation zones with a capacity for 12 gel cards.
Incubator for 24 gel cards in racks at 37 °C.
Digital interface screen for control and visualization of remaining incubation time and real temperature for each zone.

SPECIFICATIONS

Model	BGH1Q1
Old Model	BIFG-501
Temperature control range	37 °C
Temperature uniformity	± 0.5 °C
Temperature stability	± 0.2 °C
Temperature accuracy	< 0.5 °C
Heating time	25 °C~37 °C ≤ 5 min
Capacity	24 gel cards
Power	220 V 3 A 200 W
Dimension (LxWxH) mm	220x380x115 mm
Alt Name	Gel Card Incubator



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com