



FUNGAL GROWTH INCUBATOR BGH1R1 (BIFG-502)

FUNGAL GROWTH INCUBATOR BGH1R1

GEL CARD CENTRIFUGE

The gel card centrifuge holds automatic balance control and the user can set speed and centrifugation time. It can be suitable for many types of gel cards from different manufacturers.



Digital interface visor control and visualization of centrifuge speed and time remaining for processing.

Automatic balance control.

Pre-defined centrifugal parameters for optimum gel cards procession, acceleration, speed, braking and time.

Programmable audible alarms for end of centrifuge time periods.

SPECIFICATIONS

Model	BGH1R1
Old Model	BIFG-502
Speed range	800~1800 rpm
RCF. range	54~272xg
Capacity	12 gel cards
Timer	1~20 min
Noise	≤55 db(A)
Power	AC 100-240 V 85 W
Dimension (LxWxH) mm	440x340x170 mm
Weight	7 kg
Alt Name	Gel Card Centrifuge



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com