



FUNGAL GROWTH INCUBATOR BGH1R2 (BIFG-503)

FUNGAL GROWTH INCUBATOR BGH1R2

GEL CARD CENTRIFUGE

The gel card centrifuge holds automatic balance control and the user can set speed and centrifugation time. It can be suitable for many types of gel cards from different manufacturers.



Digital interface visor control and visualization of centrifuge speed and time remaining for processing.

Automatic balance control.

Pre-defined centrifugal parameters for optimum gel cards procession, acceleration, speed, braking and time.

Programmable audible alarms for end of centrifuge time periods.

SPECIFICATIONS

Model	BGH1R2
Old Model	BIFG-503
Speed range	500~1500 rpm
RCF. range	33~294xg
Capacity	24 gel cards
Timer	1~20 min
Noise	≤55 db(A)
Power	AC 100-240 V 85 W
Dimension (LxWxH) mm	530x410x180 mm
Weight	11.5 kg
Alt Name	Gel Card Centrifuge



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com