



VACUUM OVEN BGI-1103

VACUUM OVEN BGI-1103

TWO&THREE CHAMBERS VACUUM OVEN

LCD controller (with timing function)

Vacuum oven is designed especially for drying of material which is thermo-sensitive or decomposes and oxidative easily. It can be filled with inert gases, which is especially for a rapid drying of some compound material.



- Microprocessor controller.
 - Dual layer tempered glass door for clear observation.
 - Minimum heating time 50% less than traditional vacuum oven.
 - Interior chamber made from 304 stainless steel with a mirror polished finish making for easy cleaning.
 - Standard with vacuum pump.
 - Independent audible and visible temperature-limiting alarm system ensures experiments run safely (Option).
 - Each chamber can be heated, vacuumed, and injected with nitrogen independently without interfering with each other.
- Option:
- LED lamp

SPECIFICATIONS

Model	BGI-1103
Electrical Requirement	AC220V 50HZ
Controller	Controller
Power Consumption	3550W
Temperature Range	RT+10 ~ 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133Pa
Ambient Temperature	+5~40°C
Chamber Volume	216Lx2
Interior Dimension (WxDxH) mm	600x600x600/Chamber
Exterior Dimension (WxDxH) mm	880x960x1770
Shelves	3pcs /Chamber (Independent temp. control)
Independent temp. Controller	2 PCS
Chamber Material	304 Stainless steel
Packing Dimension	92*112*180
CBM	1.85
Net Weight	275.00
Gross weight	436
Alt Name	Vacuum Oven (2 chambers)

FEATURES

- Microprocessor controller.
- Dual layer tempered glass door for clear observation.
- Minimum heating time 50% less than traditional vacuum oven.
- Interior chamber made from 304 stainless steel with a mirror polished finish making for easy cleaning.
- Standard with vacuum pump.
- Independent audible and visible temperature-limiting alarm system ensures experiments run safely (Option).

- Each chamber can be heated, vacuumed, and injected with nitrogen independently without interfering with each other.

Option:

- LED lamp
- Inert gas valve
- Oil mist filter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com