



## WATER BATH BGI-4001

## WATER BATH BGI-4001

### WATER BATH

Microprocessor controller (with timing function)

Provided for precise and constant temperature and auxiliary heating in colleges industrial and mining enterprises and scientific research departments



- ABS outer shell
  - Microprocessor control with timing function
  - Digital display
  - Audible and visible alarm for over temperature
- Option:
- Intelligent programmable temperature controller

## SPECIFICATIONS

Model	BGI-4001
Electrical Requirements	220V 50Hz
Power Consumption	500W
Temperature Range	RT+5~99°C
Temperature Stability	±0.5°C
Tracking Alarm	±2°C
Chamber Volume	5L
Interior Dimension (WxDxH, mm)	280x130x150
Exterior Dimension (WxDxH, mm)	347x201x155
Timing Range	1~5999min
Remark	Two holes
Packing Dimension	52*29*30
CBM	0.05
Net Weight	7.00
Gross weight	9
Alt Name	Water Bath



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)