



WATER BATH BGI-4004

WATER BATH BGI-4004

WATER BATH

Microprocessor controller (with timing function)

Provided for precise and constant temperature and auxiliary heating in colleges industrial and mining enterprises and scientific research departments



- ABS outer shell
 - Microprocessor control with timing function
 - Digital display
 - Audible and visible alarm for over temperature
- Option:
- Intelligent programmable temperature controller

SPECIFICATIONS

Model	BGI-4004
Electrical Requirements	220V 50Hz
Power Consumption	2000W
Temperature Range	RT+5~99°C
Temperature Stability	±0.5°C
Tracking Alarm	±2°C
Chamber Volume	24.8L
Interior Dimension (WxDxH, mm)	505x330x150
Exterior Dimension (WxDxH, mm)	496x347x165
Timing Range	1~5999min
Remark	Double lines eight holes
Packing Dimension	83*42*29
CBM	0.10
Net Weight	9.00
Gross weight	13
Alt Name	Water Bath



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com