



THERMO SHAKER INCUBATOR BGI-4401

THERMO SHAKER INCUBATOR BGI-4401

THERMO-SHAKER



Microprocessor controller (with timing function)

- LED display for Current temperature and the count down time
- Metal blocks for easy replacement, cleaning and disinfecting
- Over-temperature alarm protection
- Transparent cover to reduce the experiment error

SPECIFICATIONS

Model	BGI-4401
Temperature Range	RT+5°C~100°C
Time for temp. increasing	≤15 minutes (from 20°C to 100°C)
Time for temp. decreasing	/
Temperature Accuracy	±0.5°C
Display Resolution	0.1°C
Uniformity of the block temperature	±0.5°C
Maximum Power	100W
Temperature-segment	/
Max. Time setting	5999 min
Standard	B-54X0.5ml
Packing Dimension	33*24*34
CBM	0.03
Net Weight	3.50
Gross weight	5
Alt Name	Thermo-shaker



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com