



## VERTICAL LAMINAR AIRFLOW BGI-4702

# VERTICAL LAMINAR AIRFLOW BGI-4702

## CLEAN BENCH



- Reservation on/off control for fan & UV light
- High quality, quiet fan, adjustable air volume
- Universal electrical outlet
- An additional pre-filter below work surface helps capture large particles and prolongs HEPA filter life
- Working area is equipped with floodlights and ultraviolet sterilizing lamps
- The sliding tempered glass window with self-balancing design, can be stopped any position that customer want.

## SPECIFICATIONS

Model	BGI-4702
Num. of operator	one people two side
Laminar Flow	Vertical
Cleanliness	ISO class 5, ISO-14644-1
CFU	$\leq 0.5/\text{Dish}^* \text{ h}$ ( $\phi 90\text{mm}$ petri dish)
Air Velocity	$\geq 0.3\text{m/s}$
Noise	$\leq 62\text{dB}$
Vibration	$\leq 3\mu\text{m}$
Illumination	$\geq 300\text{Lx}$
Power input	AC220V/50Hz
Peak power	250VA
External Dimensions W x D x H(mm)	870x690x520
Internal Work Area W x D x H(mm)	1010x780x1610
Packing Dimension	209*88*159
CBM	1.89
Net Weight	130.00
Gross weight	171
Alt Name	Clean Bench

## MORE INFO

We provides high-quality, safe Vertical Laminar Clean Benches. Clean benches suck air from the room or hall space, transfer the air through a HEPA filter using a fan, and then clean the bench area with filtered air.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)