



VERTICAL LAMINAR AIRFLOW BGI-4703

VERTICAL LAMINAR AIRFLOW BGI-4703

CLEAN BENCH



- Reservation on/off control for fan & UV light
- High quality, quiet fan, adjustable air volume
- Universal electrical outlet
- An additional pre-filter below work surface helps capture large particles and prolongs HEPA filter life
- Working area is equipped with floodlights and ultraviolet sterilizing lamps
- The sliding tempered glass window with self-balancing design, can be stopped any position that customer want.

SPECIFICATIONS

Model	BGI-4703
Num. of operator	two people one side
Laminar Flow	Vertical
Cleanliness	ISO class 5, ISO-14644-1
CFU	$\leq 0.5/\text{Dish}^* \text{ h}$ ($\phi 90\text{mm}$ petri dish)
Air Velocity	$\geq 0.3\text{m/s}$
Noise	$\leq 62\text{dB}$
Vibration	$\leq 4\mu\text{m}$
Illumination	$\geq 300\text{Lx}$
Power input	AC220V/50Hz
Peak power	500VA
External Dimensions W x D x H(mm)	1360x690x520
Internal Work Area W x D x H(mm)	1500x765x1610
Packing Dimension	165*91*180
CBM	2.70
Net Weight	150.00
Gross weight	215
Alt Name	Clean Bench

MORE INFO

We provides high-quality, safe Vertical Laminar Clean Benches. Clean benches suck air from the room or hall space, transfer the air through a HEPA filter using a fan, and then clean the bench area with filtered air.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com