



HORIZONTAL LAMINAR AIRFLOW BGI-4706

HORIZONTAL LAMINAR AIRFLOW BGI-4706

CLEAN BENCH



- Reservation on/off control for fan & UV light
- High quality, quiet fan, adjustable air volume
- Universal electrical outlet
- An additional pre-filter below work surface helps capture large particles and prolongs HEPA filter life
- Working area is equipped with floodlights and ultraviolet sterilizing lamps
- The sliding tempered glass window with self-balancing design, can be stopped any position that customer want.

SPECIFICATIONS

| | |
|-----------------------------------|---|
| Model | BGI-4706 |
| Num. of operator | two people one side |
| Laminar Flow | Horizontal |
| Cleanliness | ISO class 5, ISO-14644-1 |
| CFU | $\leq 0.5/\text{Dish}^* \text{ h}$ ($\phi 90\text{mm}$ petri dish) |
| Air Velocity | $\geq 0.3\text{m/s}$ |
| Noise | $\leq 62\text{dB}$ |
| Vibration | $\leq 4\mu\text{m}$ |
| Illumination | $\geq 300\text{Lx}$ |
| Power input | AC220V/50Hz |
| Peak power | 500VA |
| External Dimensions W x D x H(mm) | 1680x480x600 |
| Internal Work Area W x D x H(mm) | 1760x700x1450 |
| Packing Dimension | 191*165*86 |
| Alt Name | Clean Bench |

MORE INFO

We provides high-quality, safe Vertical Laminar Clean Benches. Clean benches suck air from the room or hall space, transfer the air through a HEPA filter using a fan, and then clean the bench area with filtered air.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com