



## STABILITY TEST CHAMBER BGI1BH6 (BTST-1606)

# STABILITY TEST CHAMBER BGI1BH6

## LARGE MEDICINE STABILITY TESTING CHAMBER



## SPECIFICATIONS

|                               |                                                                                                                                                     |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Model                         | BGI1BH6                                                                                                                                             |
| Old Model                     | BTST-1606                                                                                                                                           |
| Temperature control range     | 0~65°C                                                                                                                                              |
| Temperature Stability         | ±0.5°C                                                                                                                                              |
| Temperature Uniformity        | ±2°C ~ ±3°C                                                                                                                                         |
| Humidity Range                | 35~95%RH                                                                                                                                            |
| Humidity Stability            | ±3%RH                                                                                                                                               |
| Illumination                  | /                                                                                                                                                   |
| Illumination difference       | /                                                                                                                                                   |
| Cooling systemcooling mode    | Two sets of Imported compressor work rotationally                                                                                                   |
| Controller                    | Programmable (touch screen)                                                                                                                         |
| Sensor                        | Temp: Pt100, Humidity: capacitance sensor                                                                                                           |
| Ambient Temperature           | RT+5~30°C                                                                                                                                           |
| Electrical Requirement        | AC220V 50Hz                                                                                                                                         |
| Power Consumption             | 3500W                                                                                                                                               |
| Operating Power               | 2800W                                                                                                                                               |
| Chamber Volume                | 1000L                                                                                                                                               |
| Interior Dimension (WxDxH,mm) | 1050x590x1650                                                                                                                                       |
| External Dimension (WxDxH,mm) | 1610x890x2000                                                                                                                                       |
| Shelves                       | 4(pcs)                                                                                                                                              |
| Safety Device                 | Compressor overheating and overpressure protection, Fan overheating protection, Over temperature protection, Overload protection, Water protections |
| Remark                        | 1. Build in mini printer2. Paperless recorder (Option)3. 1000L and 1500L are double doors                                                           |
| Alt Name                      | Large Medicine Stability Testing Chamber                                                                                                            |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)