



WATER BATH BGI1BO3

WATER BATH BGI1BO3

WATER BATH

Microprocessor controller (with timing function)

Provided for precise and constant temperature and auxiliary heating in colleges industrial and mining enterprises and scientific research departments



- ABS outer shell
 - Microprocessor control with timing function
 - Digital display
 - Audible and visible alarm for over temperature
- Option:
- Intelligent programmable temperature controller

SPECIFICATIONS

Model	BGI1BO3
Electrical Requirements	220V 50Hz
Power Consumption	1500W
Temperature Range	RT+5~99°C
Temperature Stability	±0.5°C
Tracking Alarm	±2°C
Chamber Volume	24.8L
Interior Dimension (WxDxH, mm)	505x330x150
Exterior Dimension (WxDxH, mm)	496x347x165
Timing Range	1~5999min
Remark	Double lines six holes
Alt Name	Water Bath



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com