



## SHAKING WATER BATH BGI1BP3 (BBSW-3003)

# SHAKING WATER BATH BGI1BP3

## SHAKING WATER BATH

LED Microprocessor Controller (with timing function)

Widely applicable for laboratory researches on bacteria cultivation, fermentation, hybridization, chemical and biochemical reaction, enzymes and tissues research, which have a high requirement on precision of shaking speed and temperature.



## SPECIFICATIONS

Model	BGI1BP3
Old Model	BBSW-3003
Temperature Range	RT+5~100°C
Display Resolution	0.1°C
Temperature Uniformity	±1°C
Shaking Speed Range	30~180rpm
Amplitude	30mm (Standard) or 40mm (Option)
Power Consumption	1650W
Interior Dimension (WxDxH) mm	618x310x250
Exterior Dimension (WxDxH) mm	823x350x355
Alt Name	Shaking Water Bath



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)