



WATER BATH BGI1BQ6 (BBWA-3106)

WATER BATH BGI1BQ6

WATER BATH



Microprocessor controller (with timing function)

- A stamping molding stainless steel tank, easy to clean.
- LCD screen, multiple data display with timing function, easy to operate.
- Stainless steel shelves cover heater and sensor to avoid damage during using.
- Once-forming stainless steel lid.
- Cut off heater automatically in case of lack of water, meanwhile visible and audible alarms to remind users in time.
- Independent temperature-limiting alarm system.
- temperature error alarm.
- Test tube holder can be placed. (Option)

SPECIFICATIONS

Model	BGI1BQ6
Old Model	BBWA-3106
Electrical Requirement	AC220 50Hz
Power Consumption	800W
Temperature Range	RT+5~100°C
Temperature Stability	±0.2°C
Temp Alarm	0.1°C
Interior Dimension (WxDxH) mm	300x240x200
Exterior Dimension (WxDxH) mm	480x300x480
Timing Range	1~5999min
Chamber Volume	11L
Alt Name	Water Bath



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com