



MAGNETIC STIRRER BGI1BU8

MAGNETIC STIRRER BGI1BU8

MAGNETIC STIRRER

LCD Microprocessor Controller (with timing function)

The BEING Magnetic Heated Stirrers offer a large variety of sizes and temperature ranges to accommodate your laboratory needs. Our stirrers provide a user friendly PID controller offering a large LCD color screen for easy viewing. They are constructed for durability, high performance and safety.



- Large LCD screen to display more data at same time.
- External temperature sensor for liquid, temperature range from ambient temperature to 200°C .
- Free-step speed adjustment.
- Die-cast Aluminum alloy external chamber.
- Aluminum alloy working plate.
- Over-temp alarm system, auto switch off when 470°C .
- Caution indication light when plate temperature reach 50°C .

SPECIFICATIONS

Model	BGI1BU8
Stirring capacity (H ₂ O)	15
Liquid temp range	/
Working plate temp range	/
Speed range (rpm)	100~1600 (Optional 2000)
Temperature accuracy	/
Plate dimension (mm)	180x180
Exterior dimension (WxHxD) mm	210x330x130
Power (AC Hz)	220V/50
Heating power (Max. W)	-
Input power (Max. W)	50
Weight (KG)	4
Alt Name	Magnetic Stirrer





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com