



NATURAL CONVECTION OVEN BGI1C1 (BONC-101)

NATURAL CONVECTION OVEN BGI1C1

NATURAL CONVECTION OVEN



- 304 stainless steel, mirror polishing processing, easy to clean and maintain.
- Natural convection with low noise.
- PID controller with over temperature alarm and timing function ensures precise and reliable control.

Option:

- Independent over-temperature alarm system ensures experiments running safely.
- RS 485 connector can connect computer to save the data via software.

SPECIFICATIONS

Model	BGI1C1
Old Model	BONC-101
Electrical Requirement	AC220V 50HZ
Temperature Range	RT+10 ~ 200°C
Display Resolution	0.1°C
Temperature Stability	±0.5°C
Ambient Temperature	+5 ~ 35°C
Power Consumption	850W
Chamber Volume	27L
Interior Dimension (WxDxH)mm	320x300x355
Exterior Dimension (WxDxH)mm	460x520x660
Shelves	2(pcs)
Alt Name	Natural Convection Oven

APPLICATIONS

LED Microprocessor Controller (with timing function)

Provided for desiccation, torrefaction, wax-melting and sterilization in mining industry, laboratories and scientific research institutes.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com