



HIGH TEMPERATURE CIRCULATOR BGO1B1 (BCHT-101)

HIGH TEMPERATURE CIRCULATOR BGO1B1

HIGH TEMPERATURE CIRCULATOR

We offer you robust circulators for high temperature applications equipped with a circulating pump to provide excellent temperature uniformity and safe operation. Easy to maintain with LCD displays, these circulators are user friendly and give reliable results.



Microprocessor temperature controller
Easy maintenance
LCD display with touch key panel for temperature setting
Over temperature alarm system

SPECIFICATIONS

Model	BGO1B1
Old Model	BCHT-101
Capacity	5 L
Temperature Range	RT+5°C to 95°C
Temperature Accuracy	±0.2°C
Pump Flow Rate	6 L/min
Top Opening	130x130 mm
Inner Dimension (WxDxH)	240x140x150 mm
Outer Dimension (WxDxH)	390x350x520 mm
Weight	9 kg
Alt Name	High Temperature Circulator

APPLICATIONS

Enzyme reactions, Immunoassay, Incubation, Laboratory, Research, Medical, Pharmaceuticals



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com