



## HIGH TEMPERATURE CIRCULATOR BGO1B2 (BCHT-102)

## HIGH TEMPERATURE CIRCULATOR BGO1B2

### HIGH TEMPERATURE CIRCULATOR

We offer you robust circulators for high temperature applications equipped with a circulating pump to provide excellent temperature uniformity and safe operation. Easy to maintain with LCD displays, these circulators are user friendly and give reliable results.



Microprocessor temperature controller  
Easy maintenance  
LCD display with touch key panel for temperature setting  
Over temperature alarm system

## SPECIFICATIONS

Model	BGO1B2
Old Model	BCHT-102
Capacity	15 L
Temperature Range	RT+5°C to 100°C
Temperature Accuracy	±0.05°C
Pump Flow Rate	6 L/min
Top Opening	235x160 mm
Inner Dimension (WxDxH)	300x240x200 mm
Outer Dimension (WxDxH)	430x380x570 mm
Weight	12.84 kg
Alt Name	High Temperature Circulator

## APPLICATIONS

Enzyme reactions, Immunoassay, Incubation, Laboratory, Research, Medical, Pharmaceuticals



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)