



HIGH TEMPERATURE CIRCULATOR BGO1C1 (BCHT-103)

HIGH TEMPERATURE CIRCULATOR BGO1C1

HIGH TEMPERATURE CIRCULATOR

We offer you robust circulators for high temperature applications equipped with a circulating pump to provide excellent temperature uniformity and safe operation. Easy to maintain with LCD display, these circulators are user friendly and give reliable results.



Microprocessor temperature controller
Easy maintenance
LCD display with touch key panel for temperature setting
Over temperature alarm system

SPECIFICATIONS

Model	BGO1C1
Old Model	BCHT-103
Capacity	25 L
Temperature Range	RT+5°C to 100°C
Temperature Accuracy	±0.05°C
Pump Flow Rate	6 L/min
Top Opening	235x160 mm
Inner Dimension (WxDxH)	280x250x300 mm
Outer Dimension (WxDxH)	460x440x670 mm
Weight	20 kg
Alt Name	High Temperature Circulator

APPLICATIONS

Enzyme reactions, Immunoassay, Incubation, Laboratory, Research, Medical, Pharmaceuticals



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com