



## ELECTRIC HOSPITAL BED BGW-1301

# ELECTRIC HOSPITAL BED BGW-1301

MULTI-FUNCTIONAL ELECTRIC BED



## SPECIFICATIONS

Model	BGW-1301
Material	Powder coated steel bed frame, top class ABS & PP plastic
Overall dimension	2115mm(L)*1045mm(W)*477mm - 745mm
Bed board dimension	1940mm(L)*900mm(W)
Functions	5 functions with 4 electric actuators are controlled by the hand control. Backrest: 0 - 65°; Knee gatch: 0 - 40°; Trendelenburg: 0 - 12°; Rev. Trendelenburg: 0 - 12°; HI-LOW: 745mm - 477mm
Protective railings	4 damping PP protective railings with built-in angle indicators, controlled by gas springs. When protective railings fall and retract under the bed board, can achieve seamless connection with the transfer stretcher.
Net weight	139 KG
Safe load capacity	<250 KG
Castors	360° swivel castors (Dia. 150mm), with central brake.
Battery	Battery backup for power outages.
Standard accessory	IV pole
Alt Name	Multi-functional Electric Bed



## FEATURES





**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)