



STRECHER BGW1C1 (BHBD-815)

STRECHER BGW1C1

MANUAL TRANSFER STRETCHER



1. Turn the crank to adjust height. Turn clockwise, the bed board will go up. Turn anticlockwise, the bed board will go down.
2. Clear arrow symbols to instruct the user to operate.
3. Central lockable 360° swivel castors (Dia. 150mm). Retractable 5th wheel provides effortless directional movement and cornering.

SPECIFICATIONS

Model	BGW1C1
Old Model	BHBD-815
Overall Size	2013*700 mm (L*W)
Height Range (bed board to ground)	862 - 566 mm
Bed Board	1906*610 mm (L*W)
Backrest	0 - 85°
Side rails	Smooth, fast self-lowering mechanism with dampening system; connects without gap to operation table/bed
Net Weight	75 KG
Standard accessories	Mattress, I.V. Pole
Alt Name	Manual Transfer Stretcher

FEATURES

1. Turn the crank to adjust height. Turn clockwise, the bed board will go up. Turn anticlockwise, the bed board will go down.
2. Clear arrow symbols to instruct the user to operate.
3. Central lockable 360° swivel castors (Dia. 150mm). Retractable 5th wheel provides effortless directional movement and cornering.





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com