



STRECHER BGW1F1 (BHBD-819)

STRECHER BGW1F1

HYDRAULIC TRANSFER STRETCHER



SPECIFICATIONS

Model	BGW1F1
Old Model	BHBD-819
Overall dimension	2050*745 mm (L*W)
Height range (bed board to ground)	860 - 630 mm
Bed board dimension	1932*630 mm (L*W)
Backrest	0 - 85°
Knee gatch	0 - 50°
Trend/Reverse Trend	0 - 13°
Side rails	Foldable stainless steel side rails
Bed board	Radiotranslucent bed board enables X-ray imaging on stretcher
Net Weight	105 KG
Standard accessory	I.V. Pole
Optional feature	Dia. 200 mm castors optional (Height range changes to: 910 - 680 mm)
Alt Name	Hydraulic Transfer Stretcher

FEATURES





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com