



ELECTRIC HOSPITAL BED BGW1J1 (BHBD-403)

ELECTRIC HOSPITAL BED BGW1J1

MULTI-FUNCTIONAL ELECTRIC BED



SPECIFICATIONS

Model	BGW1J1
Old Model	BHBD-403
Material	Powder coated steel bed frame, top class ABS & PP plastic
Overall dimension	2160 mm (L) * 1050 mm (W) * 450 mm - 665 mm
Bed board dimension	1940 mm (L) * 900 mm (W)
Functions	Three functions controlled by hand control, with battery backup
Backrest	0 - 70°
Knee gatch	0 - 30°
HI-LOW	665 mm - 450 mm
Net weight	130 KG
Safe load capacity	< 250 KG
Castors	360° swivel castors (Dia. 150 mm) with central brake under foot board
Standard accessory	I.V. Pole
Alt Name	Multi-functional Electric Bed



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com