



MANUAL HOSPITAL BED BGW1K1 (BHBD-509)

MANUAL HOSPITAL BED BGW1K1

MANUAL BED



SPECIFICATIONS

Model	BGW1K1
Old Model	BHBD-509
Material	Powder coated steel bed frame, top class ABS & PP plastic
Overall dimension	2080 mm (L) * 980 mm (W) * 470 mm (H)
Bed board dimension	1910 mm (L) * 830 mm (W)
Functions	Two functions controlled by cranks
Backrest	0 - 75°
Knee gatch	0 - 45°
Height (bed board to ground)	470 mm
Protective railings	Collapsible 5-column aluminum protective railings, easy to grasp
Head & foot boards	PP head & foot boards, integrally blow-molded, easy to detach and clean
Castors	360° swivel castors (Dia. 125 mm) with individual brakes
Optional feature	Shoe rack under bed board (optional)
Net weight	65 KG
Standard accessory	I.V. Pole
Alt Name	Manual Bed



FEATURES





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com