



MANUAL HOSPITAL BED BGW1L1 (BHBD-510)

MANUAL HOSPITAL BED BGW1L1

3 CRANKS MANUAL BED



SPECIFICATIONS

Model	BGW1L1
Old Model	BHBD-510
Material	Powder coated steel bed frame, top class ABS & PP plastic
Overall dimension	2160 mm (L) * 970 mm (W) * 575 mm - 750 mm
Bed board dimension	1930 mm (L) * 830 mm (W)
Functions	Three functions controlled by cranks
Backrest	0 - 65°
Knee gatch	0 - 45°
HI-LOW	750 mm - 575 mm
Protective railings	Collapsible 6-column aluminum protective railings, easy to grasp
Head & foot boards	PP head & foot boards, integrally blow-molded, easy to detach and clean
Castors	360° swivel castors (Dia. 125 mm) with individual brakes
Net weight	102 KG
Standard accessory	I.V. Pole
Alt Name	3 Cranks Manual bed

FEATURES





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com