



MANUAL HOSPITAL BED BGW1M1 (BHBD-511)

MANUAL HOSPITAL BED BGW1M1

4 CRANKS MANUAL BED



SPECIFICATIONS

Model	BGW1M1
Old Model	BHBD-511
Material	Powder coated steel bed frame, top class ABS & PP plastic
Overall dimension	2120 mm (L) * 1055 mm (W) * 505 mm - 665 mm
Bed board dimension	1930 mm (L) * 900 mm (W)
Functions	Five functions activated by cranks
Backrest	0 - 60°
Knee gatch	0 - 30°
Trend/Reverse Trend	0 - 7°
HI-LOW	665 mm - 505 mm
Head & foot boards	PP head & foot boards, integrally blow-molded, easy to detach and clean
Protective railings	4 damping PP protective railings, controlled by gas springs, seamless connection with stretcher
Castors	360° swivel castors (Dia. 150 mm) with central brake under foot board
Net weight	117 KG
Standard accessory	I.V. Pole
Alt Name	4 Cranks manual bed

FEATURES





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com