



CHILDREN MEDICAL BED BGX1B1 (BHBD-104)

CHILDREN MEDICAL BED BGX1B1

MANUAL 3 FUNCTION CHILDREN MEDICAL BED

Children's medical beds are built to the most current standards while also meeting the clinician's needs which allows easy care and disinfection. The bed is equipped with Trendelenburg and reverse Trendelenburg positions. Simple and modern design for safety and protection of the child.



ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, pink, with patient name card holder, option: cartoon sticker

Full length powder coated side rail, drop type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Central locking ABS foot pedal with casters dia. 5"

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes.

SPECIFICATIONS

| | |
|-------------------------|--|
| Model | BGX1B1 |
| Old Model | BHBD-104 |
| Function | 3 |
| Back-rest lifting angle | 0°-80°±5° |
| Foot-rest lifting angle | 0°-40°±5° |
| Height adjustment | 480~720 mm (without mattress) |
| IV pole | Yes |
| Dimension | 1980Lx980Wx480~720 mm |
| Alt Name | Manual 3 Function Children Medical Bed |

APPLICATIONS

Medical unit, Nursing.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com