



ELECTRIC HOSPITAL BED BGX1H1 (BHBD-410)

ELECTRIC HOSPITAL BED BGX1H1

ELECTRIC 3 FUNCTION MEDICAL BED

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital beds allow height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.



ABS head & foot boards can be removable. With bumper wheels, 12 kg/pair, color option: blue, green with patient name card holder
2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock
Bed sheet top made from perforated steel, epoxy powder coated and baked finish
Bed framework made from mild steel ERW rectangular tubes
Central locking S.S foot pedal with casters dia. 5"
With 3 motors from Taiwan. With a battery.
Backrest & knee-rest up/down at same time
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes
Optional Accessories: Manual CPR/electric CPR.

SPECIFICATIONS

Model	BGX1H1
Old Model	BHBD-410
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	430-720 mm (without mattress)
IV pole	Yes
Dimension	2160Lx980Wx430-720 mm
Alt Name	Electric 3 Function Medical Bed

APPLICATIONS

Medical unit, ICU.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com