



ELECTRIC HOSPITAL BED BGX1J1 (BHBD-412)

ELECTRIC HOSPITAL BED BGX1J1

ELECTRIC 5 FUNCTION MEDICAL BED

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital beds allow height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.



ABS head & foot boards can be removable. 6.3 kg/pair, color options: blue with patient name card holder
ABS side rail, 4pcs, tuck away type with safe lock-out
Bed sheet top made from perforated steel, epoxy powder coated and baked finish
Bed framework made from mild steel ERW rectangular tubes
Central locking S.S foot pedal with casters dia. 5"
With 4 motors from Taiwan
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes.

SPECIFICATIONS

Model	BGX1J1
Old Model	BHBD-412
Function	5
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Trendelenburg	0°~12°
Reverse trendelenburg	0°~12°
Height adjustment	430-720 mm (without mattress)
IV pole	Yes
Dimension	2120Lx1030Wx430-720 mm
Alt Name	Electric 5 Function Medical Bed

APPLICATIONS

Medical unit, ICU.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com