



## ELECTRIC HOSPITAL BED BGX1L1 (BHBD-415)

# ELECTRIC HOSPITAL BED BGX1L1

## 5 FUNCTIONS ELECTRIC IUC BED

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital beds allow height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.



With one-key electric CPR release functions  
With manual CPR functions  
With one-key TR functions  
With one-key examination position  
With one-key chair position  
Frame: made of cold-rolled steel plate, treated by electro-coating & powder coating  
4 pcs detachable PP side rails  
Build-in control embedded inside the railing panel  
Quiet and robust electric actuators control of all reliable operation  
Nursing screen remote handset control of all electrically operated functions and movements

## SPECIFICATIONS

Model	BGX1L1
Old Model	BHBD-415
Function	5
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Trendelenburg	0°~12°
Reverse trendelenburg	0°~12°
Height adjustment	430-720 mm (without mattress)
IV pole	Yes
Dimension	2260Lx1060Wx430-720 mm
Net Weight	250 kg
Alt Name	5 Functions Electric IUC Bed

## FEATURES

With one-key electric CPR release functions  
With manual CPR functions  
With one-key TR functions  
With one-key examination position  
With one-key chair position  
Frame: made of cold-rolled steel plate, treated by electro-coating & powder coating  
4 pcs detachable PP side rails  
Build-in control embedded inside the railing panel  
Quiet and robust electric actuators control of all reliable operation  
Nursing screen remote handset control of all electrically operated functions and movements  
Pedal central-controlled caster brake system, dia. 125 mm  
Advanced medical motor system (4pcs Linak motors, 1 pcs control box, 1pcs handset)  
Accessories: 1pc telescopic SS I.V. pole; 4 pcs urine bag hooks; 4pcs I.V. holes.

## APPLICATIONS

Medical unit, ICU



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)