



ELECTRIC HOSPITAL BED BGX1N1 (BHBD-418)

ELECTRIC HOSPITAL BED BGX1N1

4 CRANK/ MANUAL 5 FUNCTION MEDICAL BED

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital beds allow height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.



ABS head & foot boards can be removable. with bumper, 8kg/pair, color option: brown, blue, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

With casters dia. 5" with individual brake

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessory: Shoe shelf, mosquito net rack, dinner table.

SPECIFICATIONS

| | |
|-------------------------|--|
| Model | BGX1N1 |
| Old Model | BHBD-418 |
| Function | 5 |
| Back-rest lifting angle | 0°-80°±5° |
| Foot-rest lifting angle | 0°-40°±5° |
| Trendelenburg | 0°~12° |
| Reverse trendelenburg | 0°~12° |
| Height adjustment | 500~780 mm (without mattress) |
| IV pole | Yes |
| Dimension | 2150Lx1060Wx500-780 mm |
| Alt Name | 4 crank/ Manual 5 function medical Bed |

APPLICATIONS

Medical unit, ICU



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com