



MANUAL HOSPITAL BED BGX1P1 (BHBD-515)

MANUAL HOSPITAL BED BGX1P1

MANUAL 2 FUNCTION MEDICAL BED

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance.



ABS head & foot boards can be removable. with bumper, 8kg/pair, color option: brown, blue, with patient name card holder
2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock.
Bed sheet top made from perforated steel, epoxy powder coated and baked finish
Bed framework made from mild steel ERW rectangular tubes
With casters dia. 5" with individual brake
Manual 2 functions adjusted by hidden type hand crank
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes
Optional Accessories: Shoe shelf, mosquito net rack, dinner table.

SPECIFICATIONS

Model	BGX1P1
Old Model	BHBD-515
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2150Lx980Wx500H mm
Alt Name	Manual 2 Function Medical Bed

APPLICATIONS

Recuperating, Nursing.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com