



## MANUAL HOSPITAL BED BGX1R1 (BHBD-517)

# MANUAL HOSPITAL BED BGX1R1

## MANUAL 2 FUNCTION MEDICAL BED

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance.



ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder  
2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock  
Bed sheet top made from perforated steel, epoxy powder coated and baked finish  
Bed framework made from mild steel ERW rectangular tubes  
Central locking ABS foot pedal with casters dia. 5"  
Manual 2 functions adjusted by hidden type hand crank  
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes  
Optional Accessories: Shoe shelf, mosquito net rack, dinner table.

## SPECIFICATIONS

Model	BGX1R1
Old Model	BHBD-517
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2160Lx980Wx500 mm
Alt Name	Manual 2 Function Medical Bed

## APPLICATIONS

Recuperating, Nursing.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)