



## MANUAL HOSPITAL BED BGY1B1 (BHBD-203)

# MANUAL HOSPITAL BED BGY1B1

## MANUAL OBSTETRIC BED

The Delivery table is a flexible designed suitable for conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.



The main bed and auxiliary bed are connected with each other, four foot wheels are equipped on auxiliary bed and brake device is spared  
Foot step hydraulic elevating, Imported Y type sealing ring and is durable  
The mattress adopts high density memory sponge and is once formed, easily cleaned and sterilized  
Rotary brake device which is easy for moving operating table  
Base seat and upright column covering are made of high quality 304 stainless steel

## SPECIFICATIONS

Model	BGY1B1
Old Model	BHBD-203
Length	1950+50 mm
Width	600+50 mm
Height Adjustment	(750~950)+50 mm
Working arealength/width	650x600 mm
Backboard upward	$\geq 50^\circ$
Backboard downward	$\geq 10^\circ$
Seat board upwards	$\geq 30^\circ$
Seat board downwards	$\geq 10^\circ$
Standard accessory	A pair of arm rods, a pair of leg support, a pair of arm rods, a pair foot pedal, a stainless steel sewage basin
Alt Name	Manual Obstetric Bed

## APPLICATIONS

Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)