



LED OPERATION LIGHT BGY1BQ1 (BOPL-603)

LED OPERATION LIGHT BGY1BQ1

LED OR SHADOWLESS LAMP

Challenging tasks require special solutions. The innovative Labstac LED operation light give you benefits of cutting-edge LED technology in a uniquely effective manner. We focus on providing LED lights that are more color-accurate, more flexible, more economical. We offer innovative light and system solutions that can be retrofitted easily and without special preparation.



Adopt imported LED cold light source, white the ultra-thin optical lens
Adopt digital means for stepless regulation of LED brightness
Adjust the diameter of light spot and color temperature
With Endoscopy model
Removable handle sheath, easy to clean and disinfect
Nano coating of heat sink has created an excellent heat dissipation effect.
Imported balance arm, which is flexible up and down

SPECIFICATIONS

Model	BGY1BQ1
Old Model	BOPL-603
Illumination	40000~160000 LUX (at 1m distance)
Color temperature	6700>TC>3000K Endo model- 10000 lux
Color rendering Index	100≥Ra≥85
Working distance	70-140 cm
Diameter of spot D50	110+35 mm
Light field diameter	10-28 cm (at 1m distance)
Heat to light ratio	<6mW/m²lx
D10	180±40 mm
Number of LED bulbs	45
Average bulb life (hours)	50000 hours
Max power consumption	90 VA
Alt Name	LED OR shadowless lamp

APPLICATIONS

Dental implant surgery, Medical unit, Surgery, Operation theatre



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com