



SINGLE-HOLE LAMP BGY1BS1 (BOPL-501)

SINGLE-HOLE LAMP BGY1BS1

SINGLE-HOLE LAMP

Single-hole lamp is suitable for urology, ENT, gynecology examination and general outpatient surgery partial auxiliary lighting. Easy to move with brakes and more suitable for cutting off power with a battery.



Streamlined, ultra-thin, enclosed head with reduced air resistance to minimise contamination

Digital controlled light intensity without grades of brightness and automatic memory for added flexibility

Optionally built in power supply systems avoiding complicated installation of wall mounted control boxes

Auxiliary light source included for laparoscopic operations

Rigid suspension arm with easy mobility and stable positioning

No infrared or ultraviolet radiation

Stepless adjusting light source

ONDAL arm, OSRAM bulb

SPECIFICATIONS

Model	BGY1BS1
Old Model	BOPL-501
Illumination	30000 LUX (at 1m distance)
Color temperature	4500K±500K
Color rendering Index	>85% Ra
Tilt left and right	180°
Tilt front and back	180°
Height adjustment distance	>400 mm
Input power	20+10% VA
Alt Name	Single-hole lamp

APPLICATIONS

Dental implant surgery, Medical unit, Surgery, Operation theatre



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com