



SINGLE-HOLE LAMP BGY1BT1 (BOPL-502)

SINGLE-HOLE LAMP BGY1BT1

SINGLE-HOLE LAMP

Single-hole lamp is suitable for urology, ENT, gynecology examination and general outpatient surgery partial auxiliary lighting. Easy to move with brakes and more suitable for cutting off power with a battery.



An imported LED cold light source is adopted
Digital means are adopted for stepless regulation of LED brightness
With electric focusing function
The removable handle sheath can be disinfected at the high temperature of 135°C.
Imported spring arm
The imported switch power supply is adopted to control the voltage, making the work voltage permanently stable
The no welding titanium alloy arm and the streamline-designed lamp cap better meet the design requirements of modern clean laminar flow operation rooms

SPECIFICATIONS

Model	BGY1BT1
Old Model	BOPL-502
Illumination	30000 LUX (at 1m distance)
Color temperature	4500K±500K
Color rendering Index	>85% Ra
Tilt left and right	180°
Tilt front and back	180°
Height adjustment distance	>400 mm
Input power	20+10% VA
Alt Name	Single-hole lamp

APPLICATIONS

Dental implant surgery, Medical unit, Surgery, Operation theatre



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com