



ELECTRIC HOSPITAL BED BGY1D1 (BHBD-205)

ELECTRIC HOSPITAL BED BGY1D1

ELECTRIC OBSTETRIC BED

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.



Imported LINAK motor various actions are controlled by remote controller
Super mute big foot wheel which is easily moved and equipped with center control breaking system.
Equipped with stainless steel sewage basin
The mattress adopts high density sponge and manufactured by high quality artificial leather without seam
With various complete accessories
Standard Accessories: A pair of leg support frame, a stainless steel sewage basin, a stand for transfusion, a remote control, a power wire

SPECIFICATIONS

Model	BGY1D1
Old Model	BHBD-205
Length	1950+50 mm
Width	800+50 mm
Height Adjustment	(550~800)+50 mm
Trendelenburg/rev-trend	≥5°/≥12°
Backboard upward	≥60°
Foot board upward	≥40°
Foot board downward	≥90°
Standard accessory	A pair of leg support, a sewage basin, a power cord, a IV pole, a remote controller
Input power	150+10%
Power Supply	220+22V, 50+1Hz
Alt Name	Electric obstetric bed

APPLICATIONS

Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com