



DELIVERY TABLE BGY1F1 (BHBD-209)

DELIVERY TABLE BGY1F1

ELECTRIC GYNECOLOGICAL TABLE

The Delivery table is a flexible design for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.



Various actions are controlled by Imported LINAK motor
Hand-held remote control and pedal remote control, independent of each other, easy to operate it
With concealed stainless steel sewage waste basin
The mattress is wrapped with imported artificial leather. High density memory sponge is formed in one piece, with no seam
Equipped with brake device which is easy for moving the table
The base is made up of high strength ABS, which is corrosion resistant and easy to clean
Leg plate can be lifted electrically
Concealed stainless steel auxiliary table
With function of auto restoration

SPECIFICATIONS

Model	BGY1F1
Old Model	BHBD-209
Table Length	1270+50 mm
Table width	550+30 mm
Height Adjustment	(700-1000)+50 mm
Trendelenburg/rev-trend	≥5°/≥12°
Backboard downward	≥40°
leg board reflex	≥30°
Working area	(350x480)+20 mm
Standard accessory	A pair of arc shaped handle, a pair of leg support, a stainless steel sewage basin, a power cord, a hand held remote controller, a foot controller
Main Voltage AC	220V+22V, 50Hz+1Hz
Alt Name	Electric gynecological table

APPLICATIONS

Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com